

Quality Control

## Certificate of Analysis

QC release date: 2020-07-30

**Product:** CK-MB FS  
**Product Code:** 11641

Cat. No.: 1 1641 99 10 021

---

**Delivery of kits:** No of kits: 30  
Kit Lot/Los: 60141085  
Kit Expiry: 2021-10

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 29015  
Expiry: 2021-10-16

Reagent 2 Lot No.: 29015  
Expiry: 2021-10-16

### Examination:

Appearance: Reagent 1: clear, colourless  
Reagent 2: clear, colourless

Odour: Reagent 1: slight smell (Thiol)  
Reagent 2: no characteristic odour

**pH:** Reagent 1 ( 25 °C ± 1 °C): 6.29  
Reagent 2 ( 25 °C ± 1 °C): 8.98  
R1+R2 ( 37 °C ± 1 °C): 6.54

**Absorbance:** Reagent 1: 0.227  
340 nm / 37 °C R1 + R2: 0.324

**Control sera:** TruLab N TruLab P  
Lot: 27547 Lot: 27548

Fresh (at 37°C) 33.1 U/L 58.5 U/L

Target values 34.7 27.8 - 41.7 63.6 50.9 - 76.3

---

**Additional remarks:** Mf-Date: 17.07.2020

This lot meets our specifications.

This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka