

Quality Control

## Certificate of Analysis

QC release date: 2020-08-17

**Product:** Creatinine PAP FS  
**Product Code:** 11759

Cat. No.: 1 1759 99 10 026

---

**Delivery of kits:** No of kits: 15  
Kit Lot/Los: 60140380  
Kit Expiry: 2021-08

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 29545  
Expiry: 2021-08-31

Reagent 2 Lot No.: 29545  
Expiry: 2021-08-31

### Examination:

Visual check Reagent 1: clear, yellow  
Reagent 2: clear, pale brown

Smell Reagent 1: like detergent  
Reagent 2: no characteristic odour

**pH:** Reagent 1 ( 25 °C ± 1 °C): 8.09  
Reagent 2 ( 25 °C ± 1 °C): 8.00

**Absorbance:** Reagent 1: 0.005  
546 nm / 37 °C R1 + R2: 0.012

**Control sera:** TruLab N TruLab P  
Lot: 25894 Lot: 25895

Test results 1.43 mg/dL 4.61 mg/dL

Target values 1.35 1.05 - 1.65 4.53 3.53 - 5.53

---

**Additional remarks:** Mf-Date: 28.06.2020

This lot meets our specifications.  
This certificate of analysis is a computer printout and has therefore not been signed  
by hand.

Sign.: Dr. G. Gorka