

Quality Control

Certificate of Analysis

QC release date: 2020-08-17

Product: Creatinine PAP FS
Product Code: 11759

Cat. No.: 1 1759 99 10 920

Delivery of kits: No of kits: 20
Kit Lot/Los: 60140382
Kit Expiry: 2021-08

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent 1 Lot No.: 29545
Expiry: 2021-08-31

Reagent 2 Lot No.: 29545
Expiry: 2021-08-31

Examination:

Visual check Reagent 1: clear, yellow
Reagent 2: clear, pale brown

Smell Reagent 1: like detergent
Reagent 2: no characteristic odour

pH: Reagent 1 (25 °C ± 1 °C): 8.09
Reagent 2 (25 °C ± 1 °C): 8.00

Absorbance: Reagent 1: 0.005
546 nm / 37 °C R1 + R2: 0.012

Control sera: TruLab N TruLab P
Lot: 25894 Lot: 25895

Test results 1.43 mg/dL 4.61 mg/dL

Target values 1.35 1.05 - 1.65 4.53 3.53 - 5.53

Additional remarks: Mf-Date: 28.06.2020

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed
by hand.

Sign.: Dr. G. Gorka