

Quality Control

Certificate of Analysis

QC release date: 2020-03-09

Product: Total Protein UC FS
Product Code: 10210

Cat. No.: 1 0210 99 10 021

Delivery of kits: No of kits: 35
Kit Lot/Los: 60138528
Kit Expiry: 2021-10

Delivery of bulk: Volume:
Processing Lot:

Results for:

Reagent Lot No.: 28549
Expiry: 2022-02-23

Examination:

Visual check: Reagent : clear, orange-red
Smell Reagent : no characteristic odour
pH: Reagent (25 °C ± 1 °C): 7.52
Absorbance: Reagent : 0.091
600 nm/ 15-25°C

Control sera: TruLab Urine 1 TruLab Urine 2
Lot: 26478 Lot: 26479
Fresh (at 37°C) 0.090 g/L 0.627 g/L
Target values 0.083 0.063 - 0.103 0.624 0.474 - 0.774

Additional remarks: Mf-Date: 25.02.2020

This lot meets our specifications.
This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka