

Quality Control

## Certificate of Analysis

QC release date: 2019-05-27

**Product:** Triglycerides FS

**Product Code:** 15710

Cat. No.: 1 5710 99 10 023

---

**Delivery of kits:** No of kits: 4  
Kit Lot/Los: 60132419  
Kit Expiry: 2021-05

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent Lot No.: 26820  
Expiry: 2021-05

### Examination:

Visual check: Reagent : clear, yellow

Smell Reagent : like phenol

**pH:** Reagent ( 25 °C ± 1 °C): 7.16

**Absorbance:** Reagent : 0.018  
546 nm/ 15-25°C

Time to reach endpoint ≤ 10 minutes at 37 °C

<b>Control sera:</b>	TruLab N Lot: 24382	TruLab P Lot: 24383
Fresh (at 37°C)	91.4 mg/dL	167 mg/dL
Target values	92.5 75.8 - 109	165 136 - 195

---

**Additional remarks:** Mf-Date: 16.05.2019

This lot meets our specifications.

This certificate of analysis is a computer printout and has therefore not been signed by hand.

Sign.: Dr. G. Gorka