

Quality Control

## Certificate of Analysis

QC release date: 2020-01-21

**Product:** NEFA FS  
**Product Code:** 15781

Cat. No.: 1 5781 99 10 935

---

**Delivery of kits:** No of kits: 4  
Kit Lot/Los: 60137393  
Kit Expiry: 2021-01

---

**Delivery of bulk:** Volume:  
Processing Lot:

---

### Results for:

Reagent 1 Lot No.: 28342  
Expiry: 2021-01-12

Reagent 2 Lot No.: 28342  
Expiry: 2021-01-12

### Examination:

Visual check Reagent 1: clear, colourless to brownish  
Reagent 2: clear, yellow to brown

Smell Reagent 1: no characteristic odour  
Reagent 2: no characteristic odour

**pH:** Reagent 1 ( 25 °C ± 1 °C): 6.98  
Reagent 2 ( 25 °C ± 1 °C): 6.97

**Absorbance:** Reagent 1: 0.003  
546 nm / 37 °C Reagent 2: 0.039  
R1+R2: 0.016

**Control sera:** TruLab L 1 TruLab L 2  
Lot: 27037 Lot: 27038

Test results 0.367 mmol/L 1.78 mmol/L

Target values 0.369 0.295 - 0.443 1.962 1.570 - 2.354

---

**Additional remarks:** Mf-Date: 13.01.2020

This lot meets our specifications.  
This certificate of analysis is a computer printout and has therefore not been signed  
by hand.

Sign.: Dr. G. Gorka